

EASY SEAT 02

Comfort cushion

Ref: TCA 09

INSTRUCTIONS FOR USE



Dear Customer!

Thank you for the confidence you have shown in purchasing Triamed products! We would like to inform you that Triamed products are made of the best performing materials, all having quality certificates and biocompatibility tests.

These devices meet the essential requirements of Government's Decision 54/2009 and have registration Certificates from the Ministry of Health.

These products satisfy the European Medical Device Directive 93/42/CEE (RD.1591/2009).

The cushion is used according to the recommendations of the specialized physician. The product must be adjusted by the orthopaedic technician!

Indications

It is used to prevent pressure sores on the pelvic area for bedridden people and to optimize patient comfort.

Description

Cushion with textile cover, filled with polystyrene beads which allow air circulation.

It reconfigures quickly, changing the pressure areas.

Dimensions: 42 cm x 42 cm x 3 cm.

Composition

Exterior: 100% cotton, interior: polystyrene beads.

Positioning

The pillow is placed under the body segment exposed to the pressure sore.



Cleaning

To be washed by hand with warm water, no bleaching, no spinning, no ironing, no chemical cleaning!

Storage

At room temperature, in a dry space with good ventilation.

The legal guarantee of conformity according to Law no. 449/2003: **2 years.**

Commercial warranty of the product: 30 days.

Date of production: see packaging label.

Validity date: no expiry date.

To validate the warranty, the seller must indicate the date of purchase. The buyer is obliged to keep and present the warranty coupon. In order to solve the problems that appeared in the warranty, the following data must be notified: Article reference, batch and date of manufacture, indicated on the packaging box of the product.

Deadline for communication of the solution: 15 days.

Please send notifications and complaints to the following email address: suportclienti@triamed.ro

WARRANTY COUPON 30 days from the date of purchase	Batch and date of manufacture
Product reference	Date of purchase
The reason of complaint	Seller signature and stamp



Manufacturer: TRIAMED SRL

România/Timiș/Timișoara, Str. Calea Buziașului nr.77, Hala 1-C1/A

Tel/Fax: 00 40 256 487388, Tel: 00 40 765 240037

office@triamed.ro, triamed@gmail.ro

www.triamed.ro