

MATERNITY

Back support

Ref. SRT 103

Indications

It is used during pregnancy giving moderate stability to the abdomen and the lumbar spine. To be used only with medical recommendation.

Description

The support has 2 components:

- the lumbar part made of elastic fabric, provided with 4 integrated plastic rods
- the frontal part made of soft, rigid fabric.

It supports the abdomen and it is adjustable so it can be used for the whole period of the pregnancy.

Hook and loop fastener on the front side.


Width: front 11 cm, back 16 cm.

For an easy application it is provided with helping strips.

The anterior rigid component provides support to the abdomen, the elastic back component stabilizes the lumbar region.


Composition

Elastic fabric: 50% polyester, 33% polyamide, 17% elastane, with a structure that allows high air permeability in both directions.




air transfer

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



air transfer



water transfer



- Available colors: black, beige

Available sizes

Universal size

For hips with circumference between 80 - 140 cm



Product details



- hook and loop closure on the front
- large adjustment surface
- straps for an easy application



- lumbar part made of elastic fabric, provided with splints for stability



- elastic fabric
- rigid textile