

CERVITEX JUNIOR

Soft cervical orthosis

Ref. SRT 125

Indications

It is used in case of cervical arthroses, torticollis, light cricks, post surgery.

Description

The collar is made of a 3 cm wide polyurethane foam (sponge) with cotton covering. It is anatomically shaped for optimal fit and increased comfort. Hook and loop fastener on the back side.

Due to the anatomical shape and to the density of the foam, the soft cervical orthosis limits the flexion and extension movements of the cervical spine. This provides support for the cervical spine, helping the healing process.

Composition

Polyurethane foam – sponge.

The fabric (covering): 92% cotton, 8% elastan, is soft and it does not cause irritation.

It allows a good ventilation, it absorbs and easily removes moisture from the body, keeping the body dry.



also available in adult version



Product details



- easy application due to the anatomical shape
- hook and loop fastening on the back side



- Available colors: black, white, grey melange

Available sizes in 2 different heights: 6 and 8 cm

Sizes	1	2	3
Neck circumference (cm)	21-24	24-27	27-30
Reference	SRT125M1H	SRT125M2H	SRT125M3H



After the Reference must be added **the height in cm: 6 or 8**

The circumference of the neck is measured