


ANKLEFIX SPORT

Mobile ankle support with metal splints

Ref. SRT 358

INDICATIONS

It is used in cricks, luxations, post consolidated fractures, inflammatory processes, after the removal of the orthopedic plaster cast.



MEDICAL



SPORT



DAYLY

DESCRIPTION

Orthosis with anatomical design, made of textile, provided with:

- 2 elastic straps crossed around the ankle for compression and fastening, which can be repositioned as needed
- 2 metal side splints one on each side
- 2 flexible metal splints, one on each side. For a perfect fit at the malleolus, the 2 splints from both sides can be switched (the metal one with the flexible one)
- the 2 elastic straps are provided with 2 pockets for the metal splints. These elastic straps can be removed or repositioned according to ankle circumference . The heel area is open.


Hook and loop fastener on the front side.

Height: 23 cm.


The orthosis is an ankle support, the 2 flexible metal splints and the 2 metal side splints ensure stability and safety in use. The crossed straps around the ankle and around the 4 splints exert additional compression, making the orthosis an efficient rigid support.

COMPOSITION

The protection pads are manufactured by thermoforming and have a shape that follows the anatomical shape of the ankle.



AIR FLOW



WATER TRANSFER

AVAILABLE COLOURS

Black



AVAILABLE SIZES

Universal size

For ankle with circumference between 20 - 30 cm.

**ONE
SIZE**

