

EASY SEAT 02

Comfort cushion

Code: TCA 09

Indications

It is used to prevent pressure sores on the pelvic area for bedridden people and to optimize patient comfort.

Description

Cushion with textile cover, filled with polystyrene beads which allow air circulation.

It reconfigures quickly, changing the pressure areas.

Dimensions: 42 cm x 42 cm x 3 cm.

Composition

Exterior: 100% cotton, interior: polystyrene beads.

Properties of the fabric

The cotton fabric is soft and does not cause irritation.

Allows good ventilation and keeps the body dry.

Absorbs and easily removes body moisture.

Extruded polystyrene beads allow free air circulation, do not absorb water vapor or moisture and have high pressure resistance. These are extremely easy to reconfigure, which leads to the reconfiguration of the pressure points, thereby helping to prevent pressure sores.



Product details



- polystyrene beads