

AKTIV LIGHT 01 JUNIOR

Mobile knee brace

Ref. SRT 333

Indications

It is used in the crick, luxation, post-fracture, inflammatory processes, and lateral knee joint instability. It is recommended both as orthopedic treatment and postoperative.

Description

It is a tubular, closed orthosis, made of textile, provided with:

- 2 spring type metal splints, one on each side
- 2 hook and loop straps for stabilization
- patellar support.

Length: varies between 22 - 25 cm, depending on the size.

The orthosis is a support for the knee joint, the 2 spring type metal splints, one for each side providing stability and safety in use.

Due to it's tubular, sleeve shape, it is not recommended to use it in case of a knee with major swelling and oedema.

Composition

Neoprene, with textile liner on both sides, 100% polyester.

Properties of the fabric

- The neoprene microporous cellular structure, with excellent shock absorption and heat storage properties, provides perfect thermal protection. Thermal effects help improve blood circulation and relaxation of the muscles.
- Due to the elastic texture of the material, neoprene products fix and protect the areas they cover, giving them flexibility. Without limiting the movements of the joints and muscles, the neoprene ensures a uniform distribution of pressure during the movement.



also available in adult version, page 53




Product details



- tubular, closed orthosis




- patellar support made of polyethylene foam



- spring-type metal splints

Available sizes

Sizes	1	2	3	4	5
Knee circumference (cm)	19-22	22-25	25-28	28-31	31-33
Height (cm)	22	22	23.5	25	25
Reference	SRT333M1	SRT333M2	SRT333M3	SRT333M4	SRT333M5



The circumference of the knee is measured at the middle of the knee cap