

TRIAFIX ARM 01

Arm support

Ref. SRT 208

Indications

It is used to support the arm and forearm.


Description

The support is made of one single textile piece that supports the arm, provided with 2 adjustable support straps. Hook and loop fastener.


The orthosis is an efficient support for the arm and forearm. Due to the two fastening straps, it does not allow the arm abduction and external rotation.

Composition

Three-layer laminated textile. The laminated textile ensemble behaves like a moisture regulator, ensuring a quick transfer of the moisture to outward.



air transfer



water transfer



Product details



- provided with adjustable support straps



- hook and loop fastener



- Available color: black

Available sizes

Sizes	1	2	3
Forearm length (cm)	22-28	28-32	32-36
Reference	SRT208M1	SRT208M2	SRT208M3

The length of the forearm is measured

